



PROJECT EPIONE

THE WEARABLE OF TOMORROW
BY GAIP IOTW 19 G1

40%

Of elderly Singaporeans
have experienced joint
pains

10%

Actually get their
injuries/pain treated

50+%

Of people are afraid of
undergoing surgeries
to cure their pain

We aim to improve
the quality of life of people,
particularly the elderly, by providing
them with guidance for simple
exercises, physiotherapeutic
movements and implement a fitness
and health monitoring system.


1

THE PROBLEM
AT HAND



2

PROTOTYPE
AND ITS
CONSTRUCTION



3

FINAL PRODUCT
AND ITS
FEATURES



4

BUSINESS
MODEL AND
EXPENSES




PROBLEM STATEMENT

To provide smart physiotherapeutic aid for the elderly.






OUR
SOLUTION

- 
1. Ensure that the exercise routine is being followed correctly in order to prevent pain and injury
 2. Monitor patient vitals and activity (steps, movements) remotely.

THE PROTOTYPE





THE NETWORK CONNECTIONS



MOBILE APPLICATION TO CLOUD

Wi-Fi

SENSOR TAG/MICROBIT TO RASPBERRY PI

Serial (wired)

RASPBERRY PI TO IOT HUB

Wi-Fi

SENSORS ON THE WEARABLE



- Temperature sensor
- Accelerometer
- Gyroscope
- Light Sensor

IOT SECURITY



DEVICE TO PI





- Encryption

PI TO HUB

- Symmetric Key Encryption

PHYSICAL SECURITY FOR HARDWARE

- Blocking the ports of the Raspberry Pi
- Smartphones have an in-built security system



THE PRODUCT





PORTABILITY

Portability can be improved tremendously by using a more compact processor (e.g: Qualcomm Snapdragon 429) instead of the Raspberry Pi



WRIST WEARABLE VARIANT

A wrist wearable variant can be introduced to cater to the requirements of different customers



ADD MORE EXERCISES

The system can be trained to detect errors in more varieties of exercises. This can be extended to yoga and aerobics as well.




INCLUSION OF A HEART RATE SENSOR

For better monitoring of vitals



BUSINESS MODEL



TYPES OF INJURIES

TRAINING INJURIES

Training errors can occur when you take on too much physical activity too quickly.

Our vital monitoring system summarizes the user's workout intensity, so that the changes can be implemented when needed.

TECHNIQUE INJURY

Improper form can lead to imbalanced loading on specific muscles and tissues, which can lead to long-term damage.

Our real time exercise monitoring helps in correcting possible faults the user could be making to prevent sustained damage.

~3

SESSIONS PER MONTH

After which the patient would have to maintain an unsupervised exercise routine.

135
SGD

COST PER SESSION

A session lasts for typically 45 mins.

CUSTOMER EXPENDITURE

Source:

- <https://www.coreconcepts.com.sg/service-fees/>
- <https://www.arthritis.org/living-with-arthritis/treatments/natural/other-therapies/what-is-physical-therapy.php>

CAPITAL EXPENSES

SENSORS

Heart rate, Temperature,
Gyroscope
 $\sim 10 + 8 + 4 = 22$ SGD

CLOUD SERVICES

~ 10 SGD/MONTH

APP DEPLOYMENT

= 34 SGD

INFRASTRUCTURE REQUIREMENTS

Smartphones(2/3/4G),
Wi-Fi connection

OPERATIONAL EXPENSES

APP MAINTENANCE

CASE STUDY




John is a 63 year old gentleman living in Singapore. Recently, he has been suffering from neck and elbow pain. Upon his family's advice, he visits a physiotherapist, who recommends 3 weekly sessions for a period of two months. Each session costs John approximately **SGD 135-150**.




CASE STUDY (CONTD.)

John decides to buy our products. Considering that the product sells at **SGD 40**, John buys it and sets up his Epione App Account. His doctor can monitor his vitals and suggest various exercises depending on how his body reacts. Now, John needs to visit his physiotherapist once a month to follow up on his reports. Each appointment costs him about **SGD 60-80**.

Monthly Cost Comparison



SNEAK PEEK



A dark blue diagonal shape, resembling a large triangle or a stylized 'L' rotated 45 degrees, occupies the right half of the image. The background is a light gray. The text 'THANK YOU' is written in a bold, teal-colored, sans-serif font within the dark blue area.

THANK YOU